


LaVue - LavueController server

LavueController server allows for communication with LaVue GUI via tango interface, e.g. with user scripts or macros



Attributes

The user can **pass to** and **read back** from the LaVue GUI the following attributes:

- **BeamCenterX**: x-coordinate of the beam center in *pixels*, e.g. used in *Angle/Q Tool*
- **BeamCenterY**: y-coordinate of the beam center in *pixels*, e.g. used in *Angle/Q Tool*
- **PixelSizeX**: x-size of the detector pixel in *micrometers*, e.g. used in *Angle/Q Tool*
- **PixelSizeY**: y-size of the detector pixel in *micrometers*, e.g. used in *Angle/Q Tool*
- **DetectorDistance**: detector distance from the sample in *mm*, e.g. used in *Angle/Q Tool*
- **Energy**: beam energy in *eV*, e.g. used in *Angle/Q Tool*
- **DetectorROIs**: **JSON** dictionary with *Regions Of Interests* ranges, e.g.
{ "pilatusrois": [[67, 131, 124, 153], [67, 69, 117, 119], [134, 129, 184, 179], [125, 72, 175, 122]] }

Moreover the user can **read** from the LaVue GUI:

- **DetectorROIsValues**: **JSON** dictionary with *Regions Of Interests* sums, e.g.
{ "pilatusrois": [8167.0, 2262.0, 478.0, 1069.0] }
- **DetectorROIsParams**: **JSON** list of image transformations performed by lavue, e.g.
["transpose", "flip-left-right", "flip-up-down"]
- **ToolResults**: **JSON** dictionary with tool results, i.e. 1d diffractogram plot or position of peaks

Finally, the user can change state of lavue by **writing** to

- **LavueState**: **JSON** dictionary with lavue configuration with parameters corresponding to command-line parameters of lavue (to display them: `lavue -h`). The currently supported commands are: `source`, `configuration`, `start`, `stop`, `imagefile`, `offset`, `rangewindow`, `dsfactor`, `dsreduction`, `filters`, `mbuffer`, `maskfile`, `maskhighvalue`, `transformation`, `scaling`, `levels`, `autofactor`, `gradient`, `viewrange`, `tool`, `toolconfig`, `tangodevice`, `doordevice`, `analysisdevice`, `log`.

e.g.

```
import tango
import json

lc = tango.DeviceProxy("p09/lavuecontroller/1")

# start the Test source
lc.LavueState = json.dumps({"source": "test", "start": True})

# stop an image source
lc.LavueState = json.dumps({"stop": True})

# set the hidra source with a server configured
lc.LavueState = json.dumps({"source": "hidra", "configuration": "haspilatus1m.desy.de"})

# start Tango Attribute image source with a tango attribute from TangoTest server
# lc.LavueState = json.dumps({"source": "tangoattr", "configuration": "sys/tg_test/1
/double_image_ro", "start": True})

# set ROI tool in the tool combobox
lc.LavueState = json.dumps({"tool": "roi"})

# set display intensity levels
lc.LavueState = json.dumps({"levels": "m20,20"})

# set level auto factor
lc.LavueState = json.dumps({"autofactor": 1})

# set auto levels
lc.LavueState = json.dumps({"autofactor": ""})

# set mask values above
lc.LavueState = json.dumps({"maskhighvalue": 30000})

# set flip-up-down transformation
lc.LavueState = json.dumps({"transformation": "flip-up-down"})

# set log intensity scaling
lc.LavueState = json.dumps({"scaling": "log"})

# set eiger image source with configuration defined by alias "eiger"
# (otherwise you have to set as configuration the whole URL string)
lc.LavueState = json.dumps({"source": "http", "configuration": "eiger", "start": True})
```

You can **read** the current LaVue state via `LavueState`. To update the `viewrange` parameter write to `LavueState` an empty JSON dictionary or `{"__update__": true}`.

Properties

The **DynamicROIsValues** device property can be set to `true` or `false`. When it is true dynamic attributes with ROIs sums are added.

They names are defined by ROIs aliases from Lavue GUI.

The **DynamicROIs** device property can be set to `true` or `false`. When it is true dynamic attributes with ROIs bounds are added.

They names are defined by ROIs aliases from Lavue GUI.

The **ROIAttributesNames** device property contains names of dynamic attributes which will be created event if the corresponding ROIs aliases in Lavue GUI are missing.